


# NeuroMag Short Protocol

To find the ideal conditions, NeuroMag™ must be tested at ratios **1  $\mu\text{L}/\mu\text{g}$  DNA**, **2  $\mu\text{L}/\mu\text{g}$** , **3  $\mu\text{L}/\mu\text{g}$**  and **3.5  $\mu\text{L}/\mu\text{g}$  DNA**. For the DNA quantity, we suggest **0.5  $\mu\text{g}$**  per well in 96-well, **1  $\mu\text{g}$**  per well in 24-well and **2  $\mu\text{g}$**  per well in 6-well.


**Seed cells to be at 70% confluent the day of transfection**

1



**Prepare 4 identical tubes of DNA**

2



**96 well plate**

**24 well plate**

**6 well plate**

0.5 $\mu\text{g}$  in 50 $\mu\text{L}$  of DMEM x 4

1 $\mu\text{g}$  in 100 $\mu\text{L}$  of DMEM x 4

2 $\mu\text{g}$  in 200 $\mu\text{L}$  of DMEM x 4

**Prepare 4 tubes of NeuroMag™ (with 4 different amounts of reagent)**

3



**96 well plate**

**24 well plate**

**6 well plate**


0.5 $\mu\text{L}$ /1 $\mu\text{L}$ /1.5 $\mu\text{L}$ /1.75 $\mu\text{L}$   
in an empty microtube

1 $\mu\text{L}$ /2 $\mu\text{L}$ /3 $\mu\text{L}$ /3.5 $\mu\text{L}$   
in an empty microtube

2 $\mu\text{L}$ /4 $\mu\text{L}$ /6 $\mu\text{L}$ /7 $\mu\text{L}$   
in an empty microtube

**Mix gently each tube of DNA (step 2) to each tube of NeuroMag™ (step 3)**

4




**Incubate 20 min at room temperature**

5




**Distribute each mix onto the cells & incubate cells 20 min on the magnetic plate**

6




**Remove the cells from the magnetic plate & incubate cells for 24 to 72h at 37°C until evaluation of transgene expression**

7



**Choose the best ratio DNA: NeuroMag™**

8



**NOTES:** (1) 24H before transfection seed the cells in 6 wells of a 96-well plate, 3 wells of a 24-well plate or 3 wells of a 6-well plate in respectively 150  $\mu$ L, 400  $\mu$ L and 2 mL of complete culture medium. (2) Allow reagents to reach RT and gently vortex them before forming complexes . (3) DMEM without serum & supplement is recommended for complexes preparation. (4) For doses less than 1 $\mu$ L, dilute your reagent with deionized water.

The complete instruction manual is accessible on OZ Biosciences website:  
[www.ozbiosciences.com](http://www.ozbiosciences.com)

#### Contact information FRANCE & EUROPE

- Technical service: +33(0) 486 948 514  
[tech@ozbiosciences.com](mailto:tech@ozbiosciences.com)
- Customer service: +33(0) 486 948 516  
[contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)

#### Contact information USA

- Technical service: +1(858) 246 7840  
[techUSA@ozbiosciences.com](mailto:techUSA@ozbiosciences.com)
- Customer service: +1(858) 246 7840  
[contactUSA@ozbiosciences.com](mailto:contactUSA@ozbiosciences.com)